

Holographic microscopy of highly focused fields via Phase Retrieval

David Maluenda^{1*}, Marcos Pérez-Aviño¹, Rosario Martínez-Herrero², and Artur Carnicer¹

¹Universitat de Barcelona (UB), Facultat de Física, Departament de Física aplicada, Martí i Franquès 1, 08028, Barcelona, Spain

²Universidad Complutense de Madrid (UCM), Departamento de Óptica, Ciudad Universitaria, 28040, Madrid, Spain

*dmaluenda@ub.edu

Abstract: Blablabla

Funding: Projects PID2019-104268GB-C21 and PID2019-104268GB-C22 by MCIN / AEI / 10.13039 / 501100011033.

```
In [1]: import os, sys
import numpy as np
import matplotlib.pyplot as plt
from pathlib import Path
import warnings
import importlib
from ipywidgets import interactive, IntSlider, Layout, HBox, VBox, TwoByTwoLayout
from IPython.display import display

warnings.filterwarnings("ignore")

home_dir = Path(os.path.expanduser('~'))
OneDrive_dir = home_dir / "OneDrive - Universitat de Barcelona"

fig_num = 0 # Figures counter
verbose = 3 # 0: Main paper figures ; 1: Also Suplementari material figures ; 2: All checks included
```

```
In [2]: IN_COLAB = 'google.colab' in sys.modules
if IN_COLAB:
    !pip install pyfftw
    !pip install -U matplotlib==3.5.2
else:
    try:
        import pyfftw
    except ImportError:
        if input("'pyfftw' module is needed. "
               "Do you want to install it via pip? "
               "[Y/n] ").lower().startswith('n'):
            print('ok... continue, so.')
        else:
            print("!pip install pyfftw")
            !pip install pyfftw
```

```
In [4]: def install_devel_module(module_name, module_git=None, module_path=None, branch=None):
    """ This function is just to import a module which might be under developing.

    There are three options:
    1. The module is already installed in the system, so let's import it.
    2. The module is not installed, but it is located in the PC, so let's install it and import it.
    3. The module is not installed, and it is not located in the PC, so let's git-clone it, install it, and imported

    :param module_name: The module that we want to import
    :param module_path: The path where to find that module. If not passed, used ./<module_name>
    :param module_git: The https route where to find that module. If not passed, used .../WavefrontEngUB/<module_name>
    :param branch: Just if you want a certain branch of this repo.
    :return: module, class, or function imported
    """

    try:
        importlib.import_module(module_name)
        print(f"'{module_name}' module already installed.")
        return
    except ImportError:
        pass

    print(f"'{module_name}' module not installed, yet.")
    module_git = "https://github.com/WavefrontEngUB/" + module_name if module_git is None else module_git
    module_path = module_name if module_path is None else str(module_path)

    do_clone = False
```

```

if 'google.colab' in sys.modules: # You are in google colab
    do_clone = True
    pip_path = module_name
    print("In Google Colab detected")
else:
    if not Path(module_path).exists():
        module_path = input(f"Type the path where you have the {module_name} module "
                             f"or left empty to git-clone it in the CWD: ")
        if module_path == "":
            do_clone = True
            pip_path = module_name
        else:
            pip_path = module_path
    else:
        pip_path = module_path
        print(f"Module {module_name} found in {module_path}")

if do_clone:
    if not Path(pip_path).exists():
        print(f"Cloning '{module_name}' module from {module_git} in {Path(pip_path).resolve()}")
        !git clone {module_git} {pip_path}
        if branch is not None:
            print(f"Checking out '{branch}' branch")
            !cd {pip_path} && git checkout {branch}

print(f"Installing '{module_name}' module from {pip_path}")
!pip install -e {"%s" % pip_path}
print(f"'{module_name}' module installed. You might need restart the kernel to import it...")

```

```

In [6]: # Getting pyHolo software
# TODO: replace this with a '!pip install -U pyHolo' or so
install_devel_module("pyHolo")

import pyHolo

```

C:\Users\dmaluenda\OneDrive - Universitat de Barcelona\Research\WorkInProgress\holographic_microscope
'pyHolo' module already installed.

```

In [7]: #Getting phase_retriever software
install_devel_module("phase_retriever")
import phase_retriever.retriever as phase_retriever

'phase_retriever' module already installed.

```

```

In [8]: pixel_size, cXv, cYv, cXh, cYh = -9999.99, 2043, 1340, 2773, 1347

```

```

In [9]: scale_ph = Path('data')/"scale.png"

w = pyHolo.misc.interactive_pixel_size(scale_ph, pixel_size, cXv, cYv, cXh, cYh, verbose=verbose, fig_num=fig_num)

```

Let's check the resolution of the system by exploring scale.png file, i.e. the effective pixel size (sampling rate).

Unckeck 'Calculate' to be more responsive when interacting with sliders.

VBox(children=(HBox(children=(VBox(children=(IntSlider(value=2043, continuous_update=False, description='Xv', ...

```

In [10]: pixel_size, fig_num = w.result
cXv, cYv, cXh, cYh = list(w.kwargs.values())[:4]

print(pixel_size, "| Box positions:", cXv, cYv, cXh, cYh) if verbose>1 else None

pixel_size.get()

```

0.037 um/pixel | Box positions: 2043 1340 2773 1347

WARNING: units are not defined, so pixel size is returned in um

```

Out[10]: 0.03650700934579439

```

Linearly polarized incident field to spherical beads of 70 nm diameter

blablabla

```

In [11]: raw_data_lin_bead40 = Path("data") / "beads_5"

lin_beads40_retriever = phase_retriever.PhaseRetriever()

lin_beads40_retriever.load_dataset(raw_data_lin_bead40, ftype='numpy')
lin_beads40_retriever.config(dim=400)
lin_beads40_retriever.config(rect=((550,550), (950,950)))
lin_beads40_retriever.config(pixel_size=pixel_size.get('lam'))
lin_beads40_retriever.config(lamb=1)
lin_beads40_retriever.config(mode='scalar')
# Lin_beads40_retriever.center_window()
lin_beads40_retriever.select_phase_origin()

# Showing the raw irradiances captured by the camera (verbose > 1)

```

```
fig_num = pyHolo.misc.plot_polarimetric_images(lin_beads40_retriever.get_images(0), "", fig_num) if verbose > 1 else fig_num
fig_num = pyHolo.misc.plot_polarimetric_images(lin_beads40_retriever.get_images(1), "", fig_num) if verbose > 1 else fig_num
```




Figure 2: beam: Polarimetric images.




Figure 3: beam: Polarimetric images.

```
In [12]: lin_beads40_retriever.compute_bandwidth(tol=4e-6)
lin_beads40_retriever.config(bandwidth=15)
if verbose > 1:
    for option in lin_beads40_retriever.options:
        print(option, lin_beads40_retriever.options[option])
```

```

pixel_size 0.07102531001127314
dim 400
rect [(550, 550), (950, 950)]
n_max 200
eps 0.01
bandwidth 15
origin (199, 184)
lamb 1
path data\beads_5
ext npy
mode scalar

```

```

In [13]: lin_beads40_retriever.retrieve()
lin_beads40_field, _ = lin_beads40_retriever.get_trans_fields()
# print(lin_beads40_field)

```

```



In [14]: # [d.keys()*d.get('scale') for d in lin_beads40_retriever.polarimetric_sets.values()]

```

```

In [15]: plt.imshow(np.abs(lin_beads40_field))
plt.colorbar()
plt.figure()
im=plt.imshow(np.angle(lin_beads40_field), vmin=-np.pi, vmax=np.pi, cmap='hsv')
cbar=plt.colorbar(im, # FIXME: The line below is not working
# cbar = fig.colorbar(im, ax=ax, cax=axis.cbar_axes[idx], orientation='vertical', shrink=0.5,
# ticks=[-np.pi, -np.pi / 2, 0, np.pi / 2,
# np.pi])
cbar.ax.set_yticklabels(
[r'$-\pi$', r'$-\pi/2$', r'$0$', r'$\pi/2$', r'$\pi$'], fontsize=20)
plt.show()

```



2. Linearly polarized incident field to spherical beads of $1\mu m$ diameter

blablabla

```
In [16]: lin_beadslum_data = Path("data") / "beads_9"

lin_beadslum_retriever = phase_retriever.PhaseRetriever()

lin_beadslum_retriever.load_dataset(lin_beadslum_data, ftype='numpy')
lin_beadslum_retriever.config(dim=200)
lin_beadslum_retriever.config(pixel_size=pixel_size.get('lam'))
lin_beadslum_retriever.config(lamb=1)
lin_beadslum_retriever.config(mode='scalar')
lin_beadslum_retriever.center_window()
lin_beadslum_retriever.config(rect=((300,300), (500,500)))
# lin_beadslum_retriever.select_phase_origin()

# Showing the raw irradiances captured by the camera (verbose > 1)
fig_num = pyHolo.misc.plot_polarimetric_images(lin_beadslum_retriever.get_images(0), "", fig_num) if verbose > 1 else fig_num
fig_num = pyHolo.misc.plot_polarimetric_images(lin_beadslum_retriever.get_images(1), "", fig_num) if verbose > 1 else fig_num
```




Figure 4: beam: Polarimetric images.






Figure 5: beam: Polarimetric images.

```
In [17]: lin_beadslum_retriever.compute_bandwidth(tol=4e-6)
lin_beadslum_retriever.config(bandwidth=25)
if verbose > 1:
    for option in lin_beadslum_retriever.options:
        print(option, lin_beadslum_retriever.options[option])

lin_beadslum_retriever.retrieve()
lin_beadslum_field, _ = lin_beadslum_retriever.get_trans_fields()
```

```
pixel_size 0.07102531001127314
dim 200
rect [(300, 300), (500, 500)]
n_max 200
eps 0.01
bandwidth 25
origin None
lamb 1
path data\beads_9
ext npy
mode scalar
```

```
In [18]: plt.imshow(np.abs(lin_beadslum_field))
plt.colorbar()
plt.figure()
im=plt.imshow(np.angle(lin_beadslum_field), vmin=-np.pi, vmax=np.pi, cmap='hsv')
cbar=plt.colorbar(im, # FIXME: The line below is not working
# cbar = fig.colorbar(im, ax=ax, cax=axis.cbar_axes[idx], orientation='vertical', shrink=0.5,
# ticks=[-np.pi, -np.pi / 2, 0, np.pi / 2,
# np.pi])
cbar.ax.set_yticklabels(
[r'$-\pi$', r'$-\pi/2$', r'$0$', r'$\pi/2$', r'$\pi$'], fontsize=20)
plt.show()
```



In []:

